

**Effective DNA & RNA decontamination in one**  
**DNA-Exitus Plus™** – Without a trace  
 A must for all genetic and forensic analysis

- Fast and safe decontamination of free nucleic acids
- Effects of the components cause a very rapid non-enzymatic, non-sequence-specific degradation of **DNA** and **RNA** molecules
- Non-toxic nor even harmful for humans
- Equipment and materials are not damaged or corroded by **DNA-Exitus Plus™**
- No aggressive acids or alkaline substances
- No toxic fumes are released
- All components are biologically degradable



Product code	Product name	Applications	Pack size
A7409,0500	DNA-Exitus Plus™ IF	DNA & RNA decontamination, indicator free version	500 mL
A7409,1000RF	DNA-Exitus Plus™ IF	refill pack	1 L
A7409,2500RF	DNA-Exitus Plus™ IF	refill pack	2.5 L
A7409,0100	DNA-Exitus Plus™ IF	sample pack	100 mL
A7089,0500	DNA-Exitus Plus™	DNA & RNA decontamination, with surface trace indicator	500 mL
A7089,1000RF	DNA-Exitus Plus™	refill pack	1 L
A7089,2500RF	DNA-Exitus Plus™	refill pack	2.5 L
A7089,0100	DNA-Exitus Plus™	sample pack	100 mL
A7153,0500	RNase-Exitus Plus™	RNase decontamination	500 mL
A7153,1000RF	RNase-Exitus Plus™	refill pack	1 L
A7153,2500RF	RNase-Exitus Plus™	refill pack	2.5 L
A7600,1000	Autoclave-Exitus Plus™	DNA & RNA degradation for autoclaving processes	6 x 1 L
A9411,0025	Exitus Plus™ Activity test	DNA-Exitus Plus activity confirmation test	25 tests

For more information visit our homepage

[http://www.itwreagents.com/download\\_file/info\\_point/IP-020/en/IP-020\\_en.pdf](http://www.itwreagents.com/download_file/info_point/IP-020/en/IP-020_en.pdf)



A2022\_EN/2022108

## How to order

### Distributors

For ordering **DNA-Exitus Plus™** and all of our other high-quality products please contact your local distributor.

**ITW Reagents** has distributors in the countries as follows:

Algeria	Egypt	Mexico	South Korea
Argentina	France	Montenegro	Spain
Armenia	Germany	Morocco	Sweden
Austria	Greece	Netherlands	Switzerland
Bangladesh	Hong Kong-China	New Zealand	Taiwan
Belarus	Hungary	Norway	Thailand
Belgium	Iceland	Pakistan	Turkey
Bosnia and Herzegovina	India	Paraguay	Ukraine
Brazil	Indonesia	Peru	United Kingdom
Bulgaria	Ireland	Poland	United States
Chile	Israel	Portugal	Uruguay
China	Italy	Romania	Vietnam
Colombia	Ivory Coast	Russia	
Croatia	Japan	Saudi Arabia	
Czech Republic	Jordan	Senegal	
Denmark	Kazakhstan	Serbia	
Ecuador	Kuwait	Slovakia	
	Lithuania	Slovenia	

For further contact details of our distributors please visit

[itwreagents.com](http://itwreagents.com)

or contact our Customer Service at

[customerservice.de@itwreagents.com](mailto:customerservice.de@itwreagents.com)

[customerservice.es@itwreagents.com](mailto:customerservice.es@itwreagents.com)

[customerservice.it@itwreagents.com](mailto:customerservice.it@itwreagents.com)

### Customers - Webshop

Authorized customers from Spain, Portugal, Austria, France and Germany can order via our webshop on [itwreagents.com](http://itwreagents.com)

